# Symphonie

Aqua System **DIGITAL**

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
Dear Symphonie Aqua User,

This operating manual is designed to provide you with the safest and best possible handling. Please read the instructions completely and attentively before initial operation and carefully store them for future accessibility. We take no liability for any claims that result from inobservance of this operating manual.

Product images may vary.

## 1. Purchased parts package

---

- 1 base plate
- 1 sensor tower incl. screw set
- 1 CCU – Central Control Unit
- 1 sensor star (a star hub with 8 arms)  
incl. LPS (Local Positioning Sensor) and sensor carrier
- 1 power supply unit
- 1 operating manual
- 1 conformity declaration
- 1 certificate incl. cloud access code
- 1 Symphonie Aqua TT cylinder
- 1 digital hose package incl. base pressure gauge  
and distribution component
- 1 set of height adjusting discs

### **Not included in delivery:**

- 1 Android Tablet or Smartphone  
(must be provided on location by customer)
- 1 control unit Symphonie Aqua VC TT System

## 2. Intended use

---

The Symphonie Aqua Digital System makes it possible to determine the surface structures of residual lower limbs in 3D, while under full weight-bearing load and optimal compression.

The Symphonie Aqua Digital System is an electronic measuring system which may only be applied exclusively in combination with the Symphonie Aqua VC control unit.

### The Symphonie Aqua Digital System consists of:

- a base plate with electro-tower
- CCU
- a sensor star
- LPS Sensor

On the floor of the base plate there are pre-milled circular openings for the TT and TF cylinders, in order to securely maintain distance to the electro-tower. Do not change the distance to the electro-tower.

The Symphonie Aqua Digital algorithm computes the precise residual limb form and its volume before providing an STL file for developing/producing a prosthetic socket by means of milling or 3D printer.

Reprocessing of the determined residual limb form is not necessary. The obtained data can be directly integrated into the production process.

## 3. Warning notice

---

Do not use the Symphonie Aqua Digital System if you have discovered that it is damaged.

Test the device before each usage to ensure that the entire system is in proper order, especially the cable connections and plugs.

Use the provided **original Symphonie Aqua Digital components** exclusively.

The use of other, unauthorized components can lead to significant damage to the Symphonie Aqua Digital System.

### The following situations can lead to substantial measurement errors:

- Mobile radio frequency communication devices (mobile telephone, radio, Bluetooth devices, WIFI)
- There must be at least a 10 meter distance between one Symphonie Aqua Digital System and an additional Symphonie Aqua Digital System.
- Make certain that the sensor star cables do not run into the measuring range of the cylinder.
- The sensor star cables should not be tread upon, wound up or bent.
- The TT cylinder must always be correctly positioned on the base plate (pre-milled, circular openings).
- Do not operate the Symphonie Aqua Digital System in the vicinity of other strong magnetic fields.

The sensor star cables should not be bent, as this can lead to a break in the cable conduit.

Do not move a TT cylinder or the sensor star (Symphonie Aqua Digital System) while measuring.

**The Symphonie Aqua Digital System must not be used under any circumstances when the patient or technician is wearing an implanted pacemaker. The magnetic**

field of the Symphonie Aqua Digital System could cause a disturbance to the pacemaker's function – a life-threatening situation!

### WARNING!

#### life-threatening danger due to risk of short circuit and/or electric shock:

- Do not disconnect cables during a measurement.
- Make sure that the Symphonie Aqua Digital System does not come in contact with liquids.
- Change the Symphonie Aqua Digital System fuses only after the electrical power supply has been disconnected.
- Do not bypass any fuses.
- Do not use any aerosol sprays in the vicinity of the Symphonie Aqua Digital System.

Ship the Symphonie Aqua Digital System only in its original packaging. **We take no liability for any damages in transit.**

For cleaning the system, use only means authorized by the Romedis GmbH. Always keep the Symphonie Aqua Digital System away from water or other liquids.

Do not place any heavy objects on the Symphonie Aqua Digital System.

Be especially careful of the cables. Ensure that no one can step on the Symphonie Aqua Digital System cables. This could lead to damaging the entire electronic assembly.

Any modifications or changes made to

the Symphonie Aqua Digital System result in an immediate cancellation of the guarantee/warranty.

At regular intervals, check that all couplings are tight and that screws on the sensor tower are tightened.

## 4. CE conformity declaration

---

Please contact the manufacturer.

## 5. Disclaimer

---

The Romedis GmbH does not guarantee any precise measurement results if the system is not operated in the form and manner prescribed in this operating manual.

The embedded computing electronics and the implementation algorithms are not single fault safe.

Operation of the Symphonie Aqua Digital System is approved only with original components included in delivery. The Romedis GmbH is not responsible for measurement results made with other components.

All components have been tested regarding electric shock (short-circuit fault).

After use, disconnect the Gardena water adapter; do not leave it connected for long terms. The Romedis GmbH is not liable for water damage.

### **Note:**

If there appears to be a problem with the device that you cannot solve – or in case you want to send the shipment back – please contact our service department beforehand.

### **Contact information**

Romedis GmbH  
Georg-Wiesböck-Ring 5a  
83115 Neubeuern  
Germany  
Phone +49 80 35 - 96 78 78 96  
Fax +49 80 35 - 96 78 78 12  
[www.romedis.com](http://www.romedis.com)  
[info@romedis.com](mailto:info@romedis.com)

## 6. Assembling the Symphonie Aqua Digital System

---

First, inspect the Symphonie Aqua Digital System for possible external damages. These must be immediately reported.

Later customer claims will be disregarded.

Connect the sensor tower cable to position 1 on the CCU. (Fig. 1)

Remove the sensor star from the protective covering and place it over the sensor carrier. (Fig. 2)

Connect the two ribbon cable connectors from the sensor star to their respective positions on the CCU (pos. 2, pos. 3). (Fig. 3, Fig.4)

Unpack the power supply unit and connect the CCU to the power supply. (Fig. 5)  
Do not turn the device on yet.

Connect the TT cylinder (using the hose package and distribution component) to the Symphonie Aqua VC control unit. The distribution component is simultaneously an inlet and drain for the water in the cylinder (connection by hose directly on cylinder). (Fig. 6)

Caution: After use, disconnect the Gardena water adapter; do not leave it connected for long terms. The Romedis GmbH is not liable for water damage.

## 7. Venting the cylinder

---

For venting the Symphonie Aqua Digital System, please see the operating instructions from the VC TF & TT Systems. Digital cylinders have vertical venting.

## 8. Hose System

---

### Attachments:

Attach the hose (1) for the water intake to the water connection. (The water pressure must have a level of at least 2 bar or 29psi.) (Fig. 7)

Attach the hose (1) -- for water intake -- to the coupling (1) on the Symphonie Aqua VC control unit (Fig. 8)

Attach the hose (2) for the water drain to the coupling (2) on the VC control unit. (Fig. 9)

Attach hoses 1 & 2 (which are connected by a distribution component) to the Symphonie Aqua Digital cylinder. (Fig. 10)

Fig. 10

Attach the hose (3) for the water drain from the VC control unit to the coupling (3) and connect this, in turn, to the water drain. (Fig. 11, Fig. 12)

Abb. 11

Finally, attach the compressed air hose. (max. 6 bar)

**During the entire process, all hoses remain connected to the system.**

## 9. Filling

### Filling the Symphonie Aqua Digital cylinder:

Fig. 13

- Place the cylinder upright and open the venting screw, in order to be sure of a complete venting.
- Remove the venting screw.
- Open the water tap.
- Use only cold water for filling the Symphonie Aqua Digital cylinder!
- Slowly start filling the cylinder, using the + button (pressure increase) with water.
- **Important: Be sure to completely fill with water, without air bubbles.**
- One knows that this is achieved when water runs out of the vent hole in this position.
- After filling, be aware of the increased weight of the cylinder (approximately 30 kg).

### After the container has been filled:

- close the venting screw again,
- place the cylinder carefully on the floor in an upright position (vestibule opening facing upward).

Protect the connections as well as the components from damage and impact.

Only use the device when completely assembled.

## 10. Solid Cup and Vector Tower

1. Now, place the SAS VC TT cylinder in the corresponding circular openings for the cylinder onto the base plate.

2. For height setting, use the included compensating plates for the Symphonie Aqua Digital System. (Fig. 14)

Fig. 14

3. Choose the solid cup according to the size of the distal end of the residual limb. (Fig. 15)

Fig. 15

4. The residual limb (stump) must not come in contact with the rim of the solid cup (danger of deformation). If a solid cup that is too large is chosen, there is too much space between the residual limb and membrane. Choose the size of the bearing cushion (distal stump end) that fits the solid cup. Recheck your choice on the patient.

5. Next, determine the correct height in the system and then measure the length from the distal end of the residual limb to the desired socket height. Calculate the difference in regard to the entire length of the height-compensating discs, always beginning with the distal end cap. (Fig. 16)

Fig. 16

(System depth – residual limb length = vector tower height)

6. Screw the distal end cap, the height-compensating discs and the solid cup (including foam cushion) together.

7. Following that, empty the water completely from the cylinder, using the reset button. (Fig. 17)

Fig. 17

8. Carefully lay the screwed together length compensation assembly in the middle of the Symphonie Aqua Digital cylinder. (Fig. 18, Fig. 19)

Fig. 18

Fig. 19

9. Now, fill the system -- up to the stump end bearing cushion -- with water (+ button). (Fig. 20)

Fig. 20

## 11. Installation of Symphonie Aqua Digital App

- Now, turn the tablet on – internet access absolutely required. (WIFI connection)
- Download the Symphonie Aqua Digital app.
- Use the access code on the certificate, in order to carry out the registration in the Symphonie Digital app.
- Turn on the CCU and connect to the tablet, per Bluetooth.
- Following a successful connection, the device is now ready for operation.

### Application note

**Please notice:** For digital measurement, the length of the residual limb must be between 50mm and 380mm!

## 12. Digitalization procedure

Integrate the patient's data (girth, MPT and body weight) to the Symphonie Aqua Digital App before putting the liner on the patient.

### General data:


Abb. 21

The app calculates the water pressure to be adjusted (see VC APP). (Fig. 21)

As soon as the recommended water pressure is reached, put the liner on the patient:

These three additional measurements must be added to the data for digital measuring: (Fig. 22)

- girth at the center of the patella (measured over the liner) distal end
- distal circumference (measured over the liner)
- distance between the patella and distal residual limb end (measured over the liner)

Fig. 23

Pull the liner until under the kneecap and lay the LPS sensor onto the patella. Secure the sensor with the liner. (Fig. 23)

Secure the 10mm-sized disc centrally on the distal end as a centering guide for the sensor star.

Remove the sensor star from the sensor carrier and lay this over the cylinder. (Fig. 24)

Fig. 24

Make sure the sensor arms are straight and fully aligned. (Fig. 25)

Fig. 25

The standing patient can now place the residual limb on the middle of the sensor star.

Pay attention to the distal centering piece.

Applying only a light load, the patient steps into the system until the end of the residual limb touches the foam surface of the vector tower. Make sure that there is a light contact with the distal bearing cushion. By sinking too deeply, the pressure distribution on the limb would be substantially affected. (Fig. 26)

Fig. 26

Make sure that the stump sinks completely into the cylinder at the desired proximal socket edge (observe linear measure of the system depth). (Fig. 27)

Fig. 27

#### Begin with the first measurement:

1. Press the “+” button on the control unit; the tank will fill with water.
2. At about 400mbar, be sure that the patient now statically bears the load completely, or as well as possible.
3. Fill the cylinder (using the + button) to the pressure determined previously by the app. During the entire impression process, pressure should be kept constant.
4. **Observe the condition of the patient. The pressure level determined by the app is only an approximate value!**
5. As soon as the calculated pressure is reached, make certain that the patient maintains full weight bearing on the residual limb while in the cylinder. During the measurement, the patient must stand calmly, using a handhold.
6. Start the measurement with the Scan command on the Symphonie Aqua Digital App.

The system confirms a measurement start with an audio signal.

The measurement takes about 60 seconds.

The end of the measurement is confirmed by a signal (4 tones).

7. For complex residual limb geometries, repeat steps 2-3 of the measurement process, as needed.
8. Following the measurement, the patient reduces the bearing load on the amputated side.
9. Using the - (minus) button, let out water until there is no longer pressure in the system. The remaining water is then pumped out of the Symphonie Aqua Digital System by using the reset button.
10. **Important: In order to avoid sinking deeper into the distal bearing cushion, the patient slowly pulls the stump out of the system, as the water pressure decreases.**
11. End the scanning process in the app. Per measurement, 64 surface points and 1 reference sensor have been determined.

Depending on the complexity of the residual limb's structure, repeat the measurement process two to three times.

#### Second measurement:

- Rotate the residual limb in a clockwise direction, approximately 10° from the main arm. Have the patient again get into the cylinder and repeat steps 1-10. (Fig. 28)

- **Please take care that the LPS sensor (under the liner) remains in a constant position during the repeated measurement.**

Fig. 28

#### Third measurement:

- Repeat points 1-10 – with a further 10° rotation (Fig. 29)

Fig. 29

**Tip: If measuring (only) twice, rotate the residual limb by 20° for the second measurement. If measuring three times, rotate the residual limb by 10° for each measurement.**

For the most precise measurement result, we recommend a 3-point measure-

ment (multi-point measurement), consequently providing you with 192 surface points from which the digital software computes the most accurate and optimal socket form.

- After completing the measurement, press the finish button on the app in order to get a 3D, virtual socket model. This procedure can take 2-3 minutes, depending on the complexity of the surface structure.
- For visual inspection, the digital app provides you with a 3D moveable socket model, which will subsequently be available (use export function) as an STL file for either milling or 3D printing.
- Send the STL file to the desired production facility.

## 13. Guarantee

---

The Romedis GmbH grants a 2 year manufacturer's guarantee for the entire system, **with the exception of the following wear and tear parts:**

- sensor star, incl. all cables
- SCU
- all components from NDI
- membrane

For these, a guarantee of 9 months is valid following initial operation of the system.

The Romedis GmbH does not accept liability for country-specific internet connections.

# IMPRESSUM | IMPRINT

---

**ROMEDIS**<sup>GMBH</sup>

Georg-Wiesböck-Ring 5a

83115 Neubeuern

Germany

Tel +49 (0) 80 35 - 96 78 78 96

Fax +49 (0) 80 35 - 96 78 78 97

info@romedis.com

www.romedis.com

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Gerichtsstand Traunstein

HRB 19576

**St. - Nr.:** 156/136/60787

## Product photography:

Fotostudio Sinnesbichler

Irmgard Sinnesbichler

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83083 Riedering

Tel +49 (0) 80 36 - 301 43 94

info@sinnesbichler.de

www.sinnesbichler.de

## Layout and print:

BRUNNER Druck + Design

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83109 Großkarolinenfeld

Tel +49 (0) 80 31 - 8 111 8

Fax +49 (0) 80 31 - 89 111 9

info@brunner-druck.de

www.brunner-druck.de

Symphonie  
Aqua System **DIGITAL**

# ROMEDIS<sup>GMBH</sup>

Georg-Wiesböck-Ring 5a  
83115 Neubeuern

Tel. +49 (0) 80 35 - 96 78 78 96

Fax +49 (0) 80 35 - 96 78 78 97

info@romedis.com

www.romedis.com

Stand: 08-19