
Symphonie

Aqua System **VC^{TT}**

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
Dear Symphonie Aqua user,

Thank you for choosing the Symphonie Aqua System VC.

This innovative impression system was developed and produced for orthopaedic technology. The quality products from the Romedis GmbH undergo extensive quality assurance and documentation, as well as close-meshed tests in order to ensure flawless production.

Max. 1 bar/14,5038 psi cylinder

Symphonie
Aqua System VC

Venting
Entlüftung

1. Safety Instructions and Preventive Measures

Please read the instructions carefully before initial operation and carefully store them for future reference.

The Symphonie Aqua System VC was developed exclusively for orthopaedic technology. It is solely intended **for creating prosthetic impressions of the lower extremities, following transtibial amputation. (For knee articulation, the Aqua pad is applied.)**

Application for other purposes has not been approved. The approved, maximal size of the lower leg stump for the Symphonie Aqua System VC is 58cm. The approved, maximal bodyweight for the Symphonie Aqua System VC is 170kg. The user is responsible for the proper application and utilisation of the device, as well as accidents and dangers which occur to people or property. This device is intended for use only by qualified, trained professionals. **Before using, check the device for damages and operational capability.** Use only spare parts or original components for the Symphonie Aqua System VC which come from the Romedis GmbH. **Do not use any caustic or acidic cleaners! A**

membrane replacement is recommended after nine months.

Products from the Romedis GmbH are constantly being further enhanced. We reserve the right to change form, technology and configuration. Therefore, no claims can be made that are derived from the descriptions and illustrations in these instructions. Please adhere to country-specific standards and regulations. These instructions provide important information regarding the correct usage, initial operation, maintenance and cleaning of the Symphonie Aqua System VC.

NOTE

- Read the instructions for use carefully and completely.
- Test the device before each usage.
- For cleaning, do not use a cleanser nor objects.
- Use only a soft cloth for cleaning.
- Lubricate the membrane excessively before each use!
- Provide a handhold or supportive device
- Do not leave patient unattended during the casting procedure.
- Before filling, open the venting screw!
- Be careful – danger of slipping!!
Provide a safe and stable standing base for the device!
- Protect the compensation container, connectors and hose.
- Do not use near electric currents or fire.
- Use only for the specified purpose.

2. Important Components

The Symphonie Aqua System VC consists of:

1. hose connector (letter 1), water intake
2. attachment for hose (letter 2), general hose mounting
3. hose connector (letter 3), water drain
4. attachment Symphonie Aqua (letter 4)
5. attachment for compressed-air hose **max. 6 bar (87,02 psi)**
6. Symphonie Aqua VC-Ready (pressure container)
7. silicone membrane
8. pressure manometer
9. stopcock
10. venting screw

1 - 4 connecting hoses


14. height adjusting discs

14.1 4 sized S (3 cm)

14.2 2 sized L (6 cm)

14.3 distal endcap (7cm)

15. cushion case

15.1 Ø 10 cm, height 8 cm

15.2 Ø 12 cm, height 8 cm

15.3 Ø 14 cm, height 8 cm

15.4 Ø 16 cm, height 8 cm


16. compensation plates

17. Aqua Pad (not included in the delivery)

- 18. Symphonie Aqua VC control unit box
- 18.1 pressure manometer/
pressure in cylinder
- 18.2 negative pressure manometer/
pressure in cylinder
Max. 1 bar (14,5038 psi)
- 18.3 – button for decompression
- 18.4 + button for increase of pressure
- 18.5 reset button, pump on

3. Labelling on the Device

Venting - Entlüftung

Only for venting after filling, according to instructions.

The serial number of the pressure container can be found on the reverse or underside of the device. The serial number of the VC unit control box can be found on the underside of the device. Please keep this information accessible in case it is needed.

4. Required Material

- When filling, use only **cold tap water without additives.**
- plaster dressing
- protective sheet
- **sufficient amount of lubricating creme for membrane**
- 1 thin nylon stocking
- bucket of plaster with cold water

5. Preparation

When delivered, all system components are fully assembled and inspected for impermeability. Before operating the system, connect the VC unit control box with the 4 hoses.

Attach the hose (4) – for the water intake to the pressure container – to the connection on the Symphonie Aqua VC control unit (Fig. 1).

Fig. 1

Attach the hose (2) – for the water drain from the pressure container – to the connection on the VC control unit and attach this to the pressure container (Fig. 2).

Fig. 2

Attach the hose (1) – for the water intake – to the VC control unit connection and attach this to the water connection (Fig. 3).

(The water pressure must have a level of at least 2 bar or 29 psi.)

Fig. 3

Attach the hose (3) – for the water drain from the VC control unit – to the connection and attach this to the water drain (Fig. 4).

Fig. 4

Finally connect the compressed air hose (max. 6 bar / 87,02 psi) (Fig. 5).

Fig. 5


During the impression process, all hoses remain connected to the system.

Lay the empty cylinder at an angle, so that the venting screw points upward (Fig. 6+8).

Fig. 6 + 8

Remove the venting screw (Fig. 6) and then open the stopcock (Fig. 7).

Open the water-tap (cold water). Slowly start filling the cylinder, using the „+ button“ (pressure increase) with water.

Use only cold water!

Lay the cylinder at a slant, so that the venting screw points upward, in order to be sure of a complete venting (Fig. 8).

Important: Be sure to completely fill with water, removing all air from the system. You will know that all air has been completely purged from the system when water runs out of the vent hole in this position (Fig. 8).

After filling, be aware of the increased weight of the cylinder (approximately 30 kg).

After the container has been filled:

- close the venting screw again (Fig. 9),

- carefully place the cylinder on the floor in an upright position (vestibule opening facing upward).

Protect the connections and components from damage and impact. Only use the device when completely assembled.

The system is now ready for use (Fig. 10).

6. Application

Safety instructions: Do not leave the patient unattended during the plaster impression process!!

Clean the membrane after every application. Observe the cleaning instructions!

Do not use any caustic or acidic cleaners!

In order to obtain a correct impression, you will need:

- cold water,
- plaster dressing,
- **lubricating creme**
(recommended: Aqua Casting Lube)
- 1 nylon stocking.

For socket systems that use a liner for mounting the prosthesis, it may be advantageous to use the corresponding liner for making the impression. If necessary, use a thin film (e.g. cling film) (Fig. 11) to protect the outer textile fibres of the liner.

For pins (Shuttle Lock) Treatment, always remove the pin (danger to the membrane).

Please make sure that the silicone membrane is in sound condition before each use (Fig. 12).

Check that the patient, as well as the cylinder, are standing securely. If needed, provide a handhold or supportive device. Check the ground surface for the danger of slipping! Use the device only on non-slip surfaces.

Before use, adjust the height of the cylinder. The top edge of the cylinder should align with the height of the sitting patient's thigh.


If height adjustment is needed, lay the provided wooden shims under the cylinder.

With the help of the VC App, determine the optimal pressure.

For this, you will need:

- the patient's weight, without the prosthesis,
- the circumference underneath the patella, at the MPT height,
- the quality of the tissue (soft - middle - hard),
- as well as the patient's K Level (K1-K4),

ranging from 1 (unathletic) to 4 (athletic) (Fig. 14).

Choose the solid cup according to the size of the distal end of the stump. (Fig. 15)

The stump must not come in contact with the rim of the solid cup (danger of deformation). If a solid cup that is too large is chosen, there is too much space between the stump and membrane. Choose the size of the distal stump end bearing cushion that fits the solid cup. Recheck your choice for the patient. Next, determine the correct height in the system and then measure the length

from the distal end of the stump to the desired socket height. Calculate the difference in regard to the entire length of the height-compensating discs, always beginning with the distal end cap. (Fig. 16)

(500mm – assumed socket length = system length to be set)

Screw the distal end cap, the height-compensating discs and the solid cup (including foam cushion) together (Fig. 17).

Empty the water completely from the cylinder, using the reset button (Fig. 18).

Carefully position the spacer assembly you have screwed together into the

cylinder (photo 19).

Now fill the system – up to the stump end

bearing cushion – with water „+ button“. To avoid membrane adhesion, lubricate the membrane above the distal bearing cushion with lubricating creme.

During the application of the plaster dressing, the patient remains seated!

In order to wrap the plaster dressing onto the liner without bubbles, we recommend pulling a thin film over the nylon stocking (Fig. 20).

Wrap the stump with wet plaster dressing, without using any pressure! Be sure to evenly cover entire limb without making cast too thick. There should be no holes or gaps. Smooth plaster thoroughly.

Tip: Use cold water to prevent plaster from hardening too soon. In order to

keep the membrane as clean as possible, pull a plastic bag over the fresh plaster before the patient inserts the limb into the cylinder.

Make sure that the patient is able to bear the load of the amputated side (while standing) for the duration of the impression.

In a standing position, the patient can now place his/her stump in the middle of the cylinder opening, so that there is light contact with the distal foam cushion.

Sinking too deeply into the cushion would considerably impair the pressure distribution on the stump.

Practical tip: For bulbous stumps, use cutting aid on cast as well so that cast can be easily removed.

Please make sure that the stump sinks completely into the cylinder at the desired proximal socket edge (observe linear measure of the system depth).

Once patient is positioned, continue to fill the cylinder with water using the „+ button“ until the patient feels a slight pressure. Once a pressure of 400mbar is achieved, as indicated on the gage, instruct the patient to bear full weight onto the limb and into the cylinder. Once patient is bearing full weight, continue to fill the cylinder (using the „+ button“ to the desired pressure as determined previously by the VC app.

During the impression process, pressure should be kept constant. Observe the condition of the patient.

The pressure level determined by the app is only an approximate value!

For safety reasons, have a **handhold or supportive device** ready. The resulting normal displacement of soft tissue takes

place during the forming of the plaster impression. The plaster hardens while the limb is in a **completely weight bearing and static state**.

During the setting time, the patient can give you valuable information regarding socket volume and bony structures of the stump. This gives the practitioner valuable information to guide any modification of the plaster cast that might be needed. (Check circumference!)

After the plaster dressing has hardened, instruct the patient to reduce the load on the amputated side.

Using the „- button“, drain water until there is no longer pressure in the system. The remaining water is then pumped out of the system by using the „reset“ button.

Important: In order to avoid sinking deeper into the distal bearing cushion, tell the patient to slowly pull the stump out of the system, as the water pressure decreases.

Next, remove the film, after which the plaster impression is removed by slowly pulling it from the stump. Empty the system completely and remove the height adjusting components. Begin the process of producing a test socket.

7. Symphonie VC-App

Symphonie VC is an orthopaedic app from the Romedis GmbH which is free of charge. It assists the orthopaedic technician in calculating the recommended pressure when providing prosthetic care, using the Symphonie Aqua System VC.

How a prosthesis fits depends primarily on the form of the socket. Control, wearing comfort and adhesion are determined by the plaster impression.

The Symphonie Aqua System makes it possible, for the first time, to determine the stump's burden within the prosthetic socket under actual conditions. Due to the resulting hydrostatic pressure, individual tissue states are revealed and can be exactly determined in the socket.


Soft tissue, bony structures or scar tissue are detected. In order to achieve a proper socket fit during the first impression, the hydrostatic pressure must be regulated for each patient when using the Symphonie Aqua System VC.

The orthopaedic technician can calculate the recommended pressure from

the VC app by means of the following parameter:

- the patient's weight in kilograms
- the circumference of the stump at the MPT (mid patella tendon) in centimeter
- the firmness of the tissue (soft, middle or hard)
- the athleticism of the patient, ranging from 1 (unathletic) to 4 (athletic).

The prevailing pressure can be adapted to the recommended pressure during the treatment, with the help of the control unit.

8. Tips for Maintenance and Care

We recommend replacing the membrane every nine months.

Clean the membrane after each usage. Follow cleaning instructions! Do not use a cleanser for cleaning the membrane. Cleanse only with water and a soft cloth. Keep the System from extreme temperatures. Do not use a high-pressure cleaner.

Do not attempt to modify the device. It is only to be used for intended purposes. Misuse can result in damage to the unit and or personal/patient injury.

For injury to persons and damages to materials occurring from improper usage, the Romedis GmbH excludes any liability.

The device is solely approved for orthopaedic technology and the production of plaster impressions of the lower extremities, following transtibial amputation.

9. Storage and Transport

During transport and storage, be careful of the connections and hoses.

Do not hit, crush or place near (or into contact with) sharp objects.

Protect the compensation tank from damage. Keep the device from impact and from turning over. Protect the membrane from caustic materials, sharp or pointed objects. Before using, check that it is in safe, operating condition.

Protect the membrane from extreme hot and cold temperatures!

The VC Unit must always be properly attached to the cylinder when starting operation.

Please contact the distributor or manufacturer concerning repairs or damages.

If you have questions regarding the Symphonie Aqua System VC, you can reach us by telephone at

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